FYNCE



FYNCE Massage gun with 4 attachments / Professional Wireless Massage gun

> 04-003A FYNCE

FYNCE would first like to thank you very much for your purchase. To reward you for your purchase you get free access to 2 of our FYNCE ebooks:

- 'Thuis sporten meer dan 50 oefeningen'
- 'Sporten Warming up & Cooling down'

Get your free ebooks via the link below!

https://www.fynce.nl/handleidingen

FYNCE recommends that you always read the manual carefully before use and that you also perform the actions according to the manual.

FYNCE is available for questions, complaints or other comments. Reach our customer service:

Mail: contact@fynce.nl

ENGLISH USER MANUAL

- The drinking bottle explained Components

- Massage gun
- 2. 4 attachments
- Carry and storage bag

The massage gun (1) is a 'gun' because it 'shoots' back and forth all the time. By placing the correct attachment on the massage gun, you can relax the muscles in your body.



The massage gun can be set in 6 positions, where 1 is the slowest position and 6 is the highest position, then the headpiece on which the attachment is mounted moves the fastest.

The massage gun comes with 4 <u>attachments</u> (2). Each attachment feels different against your body, so each attachment is intended for a different part of your body. For more information, look through the enclosed ebook 'Sporten - Warming up & Cooling down'.

Everything comes in a <u>carrying and storage case</u> (3), in this case the massage gun complete with all the supplied accessories can be easily carried and stored.

IMPORTANT PRECAUTIONS

Basic safety precautions should be followed when using this product, including:

> READ ALL INSTRUCTIONS BEFORE USE

- Before use, inspect all parts for scratches, holes, tears, worn areas, discoloration or loose stitching, and other unusual issues. If something is wrong with a part, do not use the product and contact our customer service as soon as possible.
- This product is not a toy, so it is not intended to be played with. So also keep the elements away from children, pets or any other individuals who need supervision.
- Close supervision is necessary when used by or near children.
- Before use, check that all parts are correctly positioned and attached.
- Use all parts responsibly.
- Start the exercises slowly and use smooth, controlled movements to get the parts under control.
- Use all parts only and always responsibly.
 lubiX/ FYNCE cannot be held accountable and/or responsible for the possible consequences of using this product.
- Keeps the parts from being too hot (> 35 degrees Celsius) or too cold (< -15 degrees Celsius) for a long time.
- Use the parts only in the way for which they are intended.
- Only use this product in responsible and thoughtful ways. Misuse of can lead to injuries and accidents.
- Clean all elements regularly.
- Keep this product away from water while charging.
- Do not disassemble and/or attempt to modify this product.
- Avoid strong shocks or impacts.
- Only charge this product in the manner specified in the user manual.
- If the massage gun or charger smokes, emits a foul odor, heats up or shows abnormal symptoms, there may be a risk of fire or electric shock, take the following measures as soon as possible:
 - Disconnect the charger / charging cable from the wall socket when charging is in progress.
 - Turn off the massage gun.
- · Do not use the massage gun charger for objects other than this massage gun.
- Should you ever notice any abnormalities (including skin rash or irritation), stop using the massage gun immediately and consult a physician.

FYNCE lubiX

- This product is not a medical device.
- This product is not intended for use in the freezer or microwave.
- · All elements are not waterproof, so keep them out of damp areas.
- The use of this product is not recommended for medical implants such as pacemakers.
- Avoid cleaning with chemical tools. Regular cleaning with just a damp cloth is sufficient.
- Never leave the device unattended when it is in use.

THIS PRODUCT IS FOR HOUSEHOLD USE ONLY

> PACKING CONTENTS

1x Massage gun

> Supplied in black



1x Fork shape

> Supplied in black



1x Ball shape

> Supplied in black



1x Bullet shape

> Supplied in black



1x Moon shape

> Supplied in black



1x Storage and carrying case

> Supplied in black



> TIPS FOR USE

- Maintaining the correct posture is essential while relaxing and using this massage gun. Make sure you relax your muscles.
- This massage gun can also warm up or relax your muscles, good to prepare your muscles for your workout!
- Vary your levels until you find the right one, sometimes your preference for level also varies from place to place.
- Use the attachments for the designated places, of course you can also carefully experiment with them, one attachment may feel better in a completely different place than indicated by us.
- Use the included ebooks to learn more about the massage gun and the warming up and cooling down and to get inspiration for sports exercises at home!

> CLEAN AND STORE YOUR ELEMENTS

To clean the elements, we recommend wiping with a damp cloth. Do not use soap or other cleaning agents on the straps as they may damage and/or weaken the latex

We recommend storing the massage gun and other elements in a cool, dry place away from moisture, heat and direct sunlight. Use the storage and carrying case for this, of course, because everything fits in there!

SAVE THESE INSTRUCTIONS

WWW.FYNCE.NL CONTACT@FYNCE.NL